
## Gastriode D U19157

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04/08 ALF

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
### Instruction sheet

04/08 ALF

- 1 Boss
- 2 4-mm plug for connecting anode
- 3 Anode
- 4 Grid
- 5 Boss with 4-mm plug for connecting grid
- 6 Heater filament
- 7 Cathode plate
- 8 4-mm sockets for heater filament and cathode

### 1. Safety instructions

Hot cathode tubes are thin-walled, highly evacuated glass tubes. Treat them carefully as there is a risk of implosion.

- Do not subject the tube to mechanical stresses.
- Do not subject the connection leads to any tension.

If voltage or current is too high or the cathode is at the wrong temperature, it can lead to the tube becoming destroyed.

- Do not exceed the stated operating parameters.

When the tube is in operation, the terminals of the tube may be at high voltages with which it is dangerous to come into contact.

- Only change circuits with power supply equipment switched off.
- Only exchange tubes with power supply equipment switched off.

When the tube is in operation, the stock of the tube may get hot.

- Allow the tube to cool before dismantling.

The tube may only be used with tube holder D (U19100).

The compliance with the EC directive on electromagnetic compatibility is only guaranteed when using the recommended power supplies.

### 2. Description

The gas triode allows recording of the  $I_A - U_A$  characteristic of a thyatron, observation of excited and non-excited discharge and of the discontinuous transfer of energy by helium atoms in inelastic collision with free electrons.

The gas triode is a helium-filled tube with a pure tungsten heater filament (cathode) and a round metal plate (anode) with a wire grid between them, all inside a clear glass bulb. The cathode, anode and grid are all aligned parallel to one another. This planar configuration corresponds to the conventional symbol for a triode. A circular metal plate attached as a backing to the filament ensures that the electric field between the anode and cathode is uniform.

### 3. Technical data

Gas filling:	Helium
Filament voltage:	$\leq 7.5$ V AC/DC
Anode voltage:	max. 500 V DC max.
Anode current:	10 mA typ. at $U_a = 300$ V
Grid voltage:	max. 30 V
Glass bulb:	130 mm diam. approx
Length of tube:	260 mm approx.

### 4. Operation

To perform experiments using the gas triode the following equipment is also required:

1 Tube holder D	U19100
1 DC Power supply 0 – 500 V	U33000-115
or	
1 DC Power supply 0 – 500 V	U33000-230
1 Analogue multimeter AM50	U17450

#### 4.1 Setting up the tube in the tube holder

- The tube should not be mounted or removed unless all power supplies are disconnected.
- Push the jaw clamp sliders on the stanchion of the tube holder right back so that the jaws open.
- Push the bosses of the tube into the jaws.
- Push the jaw clamps forward on the stanchions to secure the tube within the jaws.

#### 4.2 Removing the tube from the tube holder

- To remove the tube, push the jaw clamps right back again and take the tube out of the jaws.

### 5. Example experiments

#### 5.1 Non-self-sustaining discharge

5.1.1 The transport of charge as compared to a vacuum triode by recording the  $I_a - U_a$  characteristics

Assuming that the number of electrons emitted by the hot cathode is identical for the same filament voltage, a comparison of the  $I_a - U_a$  characteristics to deduce the presence of the negative charge carriers that are produced when gas discharge arises as a result of ionisation due to collisions between electrons and helium atoms is possible.

- Set up the circuit as in fig. 1 and dim the lights in the room to view the faint light.
- Apply a voltage of 6 V to the heater and wait about 1 minute for the heater temperature to achieve thermal stability.

- Increase the anode voltage  $U_a$  from 0 V to 50 V in 10 V steps, then from 50 V to 250 V in 50 V steps. For each step measure the corresponding anode current  $I_a$ .

At about 25 V the anode current  $I_a$  increases considerably in the gas triode. This increase is accompanied by the appearance of a blue luminescence. It is apparent that there are many more charge carriers transporting charge than in the vacuum triode.

#### 5.1.2 Evidence of positive charge carriers

To determine the polarity of the charges detected at the grid, a sensitive calibrated measuring amplifier is used, e.g. U8531401.

- Set up the equipment as in fig. 2. Set the measuring amplifier to the range  $10^{-7}$  to  $10^{-9}$  A.
- Gradually increase the anode voltage  $U_a$  from 0 to 30 V.

The first charge carriers that reach the grid are negatively charged. At about 25 V the gas starts to discharge and the number of positive helium atoms resulting from ionisation by colliding electrons increases enormously.

#### 5.1.3 Ionisation by UV light

In this experiment gas discharge does not result from ionising collisions but from energetic light quanta (UV light) with the cathode cold. A mercury high-pressure lamp is required.

- Set up the experiment as in fig. 3. Dim the lights in the room.
- Before the start of the experiment illuminate the space between the electrodes with the UV lamp (mercury high-pressure lamp) for 1 minute.
- Gradually increase the anode voltage  $U_a$  and determine the striking potential  $U_s$  for the gas discharge (step change in current).

If the striking voltage  $U_s$  is higher than the extinguishing voltage  $U_e$  (see 5.2) then the gas discharge takes place without lighting.

#### 5.2 Determining the striking voltage and extinguishing voltage with self-sustaining discharge

In this experiment gas discharge does not result from ionising collisions but from energetic light quanta (UV light) with the cathode cold.

- Set up the experiment as in fig. 3. Dim the lights in the room.
- Gradually increase the anode voltage  $U_a$  and record the 'striking' voltage  $U_s$  for gas discharge.

- Reduce the anode voltage again until the discharge ceases. Record this extinguishing voltage  $U_E$ , too.

Another experiment can show how an excited discharge can continue without the exciting trigger.

- To this end set the anode voltage to about 10 V more than the extinguishing voltage.
- Irradiate the gas-filled space between the electrodes with UV light.
- Once the anode voltage is applied, turn off the mercury lamp.

The discharge that was excited (by UV light) still continues although the exciting trigger has been removed.

### 5.3 Inelastic collisions between electrons and helium atoms

The electrons emitted from the cathode are accelerated into the space between the cathode and the grid (held at earth potential) and, having passed through the mesh of the grid, only arrive at the anode (at -6 V) if they possess sufficient kinetic energy to overcome the field acting against them. Electrons collide with helium atoms, either elastically or, at certain kinetic energies, inelastically whereby they give up some of their energy. Electrons that overcome the field cause a current to be measured between the anode and ground.

- Set up the experiment as in fig. 4. Dim the lights in the room.
- For a reverse polarity  $U_R$  of 6 V gradually raise the accelerating potential  $U_A$  from 0 V to 70 V and observe the anode current  $I_A$ .

Up to an accelerating potential of about 24 V, the current increases (elastic collision) but then it drops suddenly. Electrons are now colliding inelastically with the He atoms and thus losing kinetic energy. This causes the helium atoms to become excited and they start to emit quanta of light.

As the accelerating potential  $U_A$  is further increased the current  $I_A$  increases once again but after another 20 V or so it drops again.

A plot of the anode current should exhibit two clear maxima. If this is not perceptible, the filament voltage should be lowered somewhat.

## 5.4 Characteristics of a thyratron

### 5.4.1 Recording the $I_A - U_A$ characteristic of a thyratron

To investigate the function of the grid, determine the anode current  $I_A$  as a function of the grid voltage  $U_G$  at constant anode voltage  $U_A$ .

- Set up the experiment as in fig. 5a. Dim the lights in the room and apply the filament voltage.
- Set the grid voltage  $U_G$  to -20 V.

- Set the anode voltage  $U_A$  to 100 V.
- Increase the grid voltage  $U_G$  to 0 V in steps of 2 V then go back to -20 V. Plot the change in the anode current  $I_A$ .

At about 10 V the gas discharge 'strikes'. Further changes to the grid have no effect on the anode current (see fig 5b) because all the positive charge carriers produced by the discharge accumulate around the grid (where the potential is lowest) and cannot pass on their charge.

### 5.4.2 Striking and extinguishing voltages as a function of the grid voltage

- Set up the experiment as in fig. 5a. Dim the lights in the room.
- Set various grid voltages  $U_G$  (-20 V, -15 V, -10 V etc.) and determine the relevant striking and extinguishing voltages by raising and lowering the anode potential  $U_A$ .

### 5.4.3 Thyratron rectifier

To use the gas triode as a rectifier it needs to be set up with an alternating anode potential from an isolating transformer. An oscilloscope is used to measure the voltage drop across a resistor that is directly proportional to the changes in the anode current.

- Set up the experiment as in fig. 6. Dim the lights in the room.
- Apply voltage to the filament.
- Apply the anode potential with a grid voltage  $U_G$  of 0 V.
- Reduce the grid voltage down to -20 V and observe the change in the anode current on the oscilloscope

### 5.4.4 Thyratron oscillator

To generate oscillations of various frequencies, a capacitor connected in parallel with the anode is charged via a resistor. When the 'striking' voltage  $U_S$  is attained, the capacitor quickly discharges. When the capacitor voltage drops below the extinguishing voltage conduction ceases and the capacitor charges back up again. The frequency of the oscillation  $f$  is not only dependent on the time constant RC but on the negative potential  $U_G$  at the grid, which affects the 'striking' voltage  $U_S$ . This means that changing the grid voltage changes the frequency of oscillation  $f$ .

- Set up the experiment as in fig. 7. Dim the lights in the room.
- Apply a voltage of -20 V to the grid.
- Perform the experiment with anode potentials between 100 V and 150 V.
- In each experiment, gradually increase the grid voltage until the gas discharge begins.

This 'striking' process is characterised by brief luminescence and the simultaneous reduction of the capacitor voltage.

- By counting the instances of luminescence and measuring the time with a stopwatch, deter-

mine the frequency of oscillation for various anode and grid potentials.

As the striking voltage decreases (i.e. increasing grid voltage) the frequency gets larger.


Fig. 1 Recording the  $I_A - U_A$  characteristic


Fig. 2 Evidence for positive charge carriers


Fig.3 Ionisation by UV light


Fig. 4 Inelastic collisions between electrons and helium atoms


Fig. 5a Recording the  $I_A - U_A$  characteristic of a thyatron


Fig. 5b


Fig. 6 Thyatron rectifier



Fig. 7 Thyatron oscillator

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